


M227A/B**RotaPill Pilling & Snagging Tester**

The RotaPill Pilling and Snagging Tester offers numerous options for testing fabrics. From the standard ICI Pilling and Snagging Box to M&S drums to the Octagonal Snagging Drum, the RotaPill will make your testing reliable and easy. The color touchscreen controller has all the standard tests preloaded so tests can be started with just a couple of touches. Magnetic backed cork liners make cork changes easy without the need to remove the box from the drive.

- Two sizes available – 2 or 4 stations
- Operates ICI Pilling/Snagging Box, Octagonal Snagging Drum, & M&S Pilling & Snagging Drum
- Drives can accommodate different boxes and drums simultaneously
- Snagging points can be mounted in the pilling boxes – no need for different boxes
- Touch screen controller has standard tests preloaded
- Count down timer lets operator know when test will finish
- Cork liners held by magnets for easy changing – no need to remove boxes
- One piece, molded ICI boxes

STANDARDS:

| | | | | | | | |
|----------------|----------|---------|----------------|-------------|------------|----------------|----------|
| adidas TM 4.08 | BS 5811 | BS 8479 | EN ISO 12945-1 | GB/T 4802.3 | ICI TM 444 | Woolmark TM152 | M&S P18A |
| M&S P18B | M&S P21A | Next 19 | | | | | |

M227A SIZE *Included 2 Pilling Boxes*

| | |
|------------|---|
| WEIGHT | 80 kg / 176 lbs |
| DIMENSIONS | 830 mm (w) x 460 mm (d) x 720 mm (h) / 33 in (2) x 18.5 in (d) x 28.5 (h) |

M227B SIZE *Included 4 Pilling Boxes*

| | |
|------------|--|
| WEIGHT | 100 kg / 220 lbs |
| DIMENSIONS | 830 mm (w) x 460 mm (d) x 935 mm (h) / 33 in (w) x 18.5 in (d) x 37 in (h) |

ORDERING INFORMATION**INSTRUMENT**

| | |
|-------|-----------------------------------|
| M227A | RotaPill Drive System - 2 station |
| M227B | RotaPill Drive System - 4 station |

ACCESSORIES

| | |
|--------|--|
| 107700 | ICI Pilling & Snagging Box with cork lined plates and four tubes |
| 101031 | EasyMount Sample Mounting Jig |
| 101035 | Standard Sample Mounting Jig |
| 101034 | Octagonal Snagging Drum for BS 8479 includes four sample tubes with locking rings (requires 200947 felt for covering sample tubes) |
| 106873 | M&S Pilling & Snagging Drum with liner, step, and 8 sample tubes with locking rings |
| 106883 | M&S Sample Template for M&S Pilling Test |
| 106891 | M&S Sample Template for M&S Snagging Test |
| 106897 | Glass Bean Bags (set of 4) for M&S Drum |
| 106894 | Pinned Snagging Bars (set of 3) to convert one M&S Drum for Snagging |

ASSESSMENT

| | |
|---------|--|
| M227G | PillGrade Automatic Pilling Grading System |
| M227C | Pillscope Assessment Viewer |
| M227PAV | Universal Pilling Assessment Viewer |
| 200999 | M&S Holoscope Assessment Viewer |
| 107608 | Pilling Photographs (set of 15) |

CONSUMABLES

| | |
|--------|---|
| 101026 | Sample Tubes (pack of 4) – for ICI Pilling Box |
| 101027 | Unmounted Cork Linings (pack of 6) – for ICI Pilling Box |
| 107701 | Snagging Points (set of 6) able to be mounted in 107700 Pilling Box |
| 106841 | M&S Sample Tubes & Locking Rings (set of 8) |
| 106885 | M&S Liner & Step for M&S Pilling & Snagging Drum |
| 200947 | Felt for BS 8478 Octagonal Snagging Drum Test (1 m x 1.5 m) |